

# Vivera Monobloc M301

**STANDARD FEATURES**

- Battery-powered hydraulic lift system
- Integrated battery with battery management system
- Automatic battery charging in fastener system
- Manual back-up release handle
- Battery level indicator
- Settable load height
- High visibility powder-coated frame
- Integrated bumper system
- Retractable head section
- One-hand release, folddown siderails
- One-hand release, infinite positioning, pneumatically assisted backrest and legrest
- Rubber wheels for better shock absorption
- Steer lock system enables wheels either to swivel or to lock into straight-line position
- One four-point shoulder restraint and two lap belts
- Anatomic mattress and pillow
- Foot-end wheel locks

**OPTIONAL FEATURES**

- Mattress with integrated children restraint
- Automatic restraints
- Head extension
- Siderails extending patient surface
- One-hand release, infinite positioning headrest (not compatible with head extension)
- One-hand release, infinite positioning kneerest



**SPECIFICATION**

Type	<b>M301</b>
Length	
Standard	<b>200 cm</b>
Minimum	<b>160 cm</b>
Width	<b>60 cm</b>
Height	
Highest position	<b>133 cm</b>
Lowest position	<b>38 cm</b>
Weight <sup>1</sup>	<b>75 kg</b>
Wheels	
Diameter	<b>15 cm</b>
Width	<b>5 cm</b>
Articulation	
Backrest	<b>0° - 90°</b>
Shock position	<b>0° - 30°</b>
Headrest <sup>2</sup>	<b>-30° - +35°</b>
Kneerest <sup>2</sup>	<b>0° - 50°</b>
Head extension <sup>2</sup>	<b>+29 cm</b>
Extended patient surface <sup>2</sup>	<b>94 cm</b>
Battery capacity	<b>5 Ah</b>
Maximum load capacity	<b>320 kg</b>
Recommended loading height	<b>&lt; 95 cm</b>
Fastener system	<b>F401</b>

<sup>1</sup> Stretcher is weighed with battery pack, without mattress and restraints.

<sup>2</sup> Optional.

Meets dynamic crash standards EN-1789 with F401 fastening system.

Specifications are rounded to the nearest whole number.

MEDIROL reserves the right to change specifications without notice.