

Clinic Extero

N114-P113

STANDARD FEATURES

- Remote steer lock system enables wheels either to swivel or to lock into straight-line position
- Seven height positions
- Automatic locking of legs when lifting
- One-person operation
- Roll-in loading with two-stage load release mechanism
- Foot-end wheel locks
- Reflective labeling
- Oversized wheels maximize mobility on all types of terrain
- Extendable one-hand release lift handles
- Light-weight aluminum construction
- Pressure-washable
- One-hand release, fold-down side rails
- One-hand release, infinite positioning, lift assisted, pneumatic backrest and legrest

OPTIONAL FEATURES

- Safety Lock Device helps eliminate accidental patient drops when unloading
- Additional lock for locking the legs in the lowest position
- Side extendable lift handles
- One-hand release, infinite positioning, lift assisted, pneumatic headrest
- One-hand release, infinite positioning, lift assisted, pneumatic kneerest

¹ Stretcher is weighed without mattress and restraints.

² Meets dynamic crash standards EN-1789 with F102 fastening system.

³ Optional.

Specifications are rounded to the nearest whole number.

MEDIROL reserves the right to change specifications without notice.



TRANSPORTER SPECIFICATION

Type	P113
Length	194 cm
Width	56 cm
Loading height	59 cm or 70 cm
Height (loading height 59/70 cm)	
Position 1	89/100 cm
Position 2	84/94 cm
Position 3	78/87 cm
Position 4	72/80 cm
Position 5	66/73 cm
Position 6	60/64 cm
Position 7	32/32 cm
Weight	28 kg
Wheels	
Diameter	20 cm
Width	5 cm
Maximum load capacity	275 kg
Fastener system ²	F102 Rolfix

STRETCHER SPECIFICATION

Type	N114
Length	195 cm
Width	56 cm
Height	12 cm
Weight ¹	20 kg
Maximum load capacity	250 kg
Handle extension	27 cm
Articulation	
Backrest	0° - 90°
Legrest	0° - 30°
Headrest ³	-30° - +35°
Kneerest ³	0° - 50°