

Vivera Clinic N114-P113

STANDARD FEATURES

- Battery-powered hydraulic lift system
- Integrated battery with battery management system
- Automatic battery charging in fastener system
- Manual back-up release handle
- Battery level indicator
- Settable load height
- High visibility powder-coated frame
- Integrated bumper system
- One-hand release, fold down siderails
- One-hand release, infinite positioning, pneumatically assisted backrest and legrest
- Rubber wheels for better shock absorption
- Steer lock system enables wheels either to swivel or to lock into straight-line position
- One four-point shoulder restraint and two lap belts
- Anatomic mattress and pillow
- Foot-end wheel locks
- External battery charging using mains power supply

OPTIONAL FEATURES

- Mattress with integrated children restraint
- Automatic restraints
- Headextension
- Siderails extending patient surface
- One-hand release, infinite positioning headrest (not compatible with head extension)
- One-hand release, infinite positioning kneerest

¹ Transporter is weighed with battery pack.

² Meets dynamic crash standards EN-1789 with F401 fastening system.

³ Stretcher is weighed without mattress and restraints.

⁴ Optional.

Specifications are rounded to the nearest whole number.

MEDIROL reserves the right to change specifications without notice.



TRANSPORTER SPECIFICATION

Type	P401
Length	200 cm
Width	60 cm
Height of seat section	
Highest position	142 cm
Lowest position	47 cm
Weight ¹	71 kg
Wheels	
Diameter	15 cm
Width	5 cm
Battery capacity	5 Ah
Maximum load capacity	300 kg
Recommended loading height	< 95 cm
Fastener system ²	F401

STRETCHER SPECIFICATION

Type	N114
Length	195 cm
Width	56 cm
Height	12 cm
Weight ³	20 kg
Maximum load capacity	250 kg
Handle extension	27 cm
Head extension ⁴	+ 29 cm
Extended patient surface ⁴	94 cm
Articulation	
Backrest	0° - 90°
Shock position	0° - 30°
Headrest ⁴	-30° - +35°
Kneerest ⁴	0° - 50°