


Vivera Fastener F401

MEDIROL

With patented technology and automated features, the Vivera Fastener loading and fastening system increases operator and patient safety when loading and unloading the powered stretchers.

Key features and benefits

- Fully mechanical loading and fastening system
- Integrated stretcher charging
- Guide system navigates stretcher during loading and unloading for maximum safety and performance
- Designed for easy installation to the ambulance floor
- Conforms to EN 1789



400 kg

MAXIMUM LOAD CAPACITY



66 kg

WEIGHT



Up to 95 cm

LOADING HEIGHT



Specifications

Type	F401
Length	200 cm
Length with extended arms	311 cm
Width	57 cm
Height	26 cm
Weight	66 kg
Maximum load capacity	400 kg
Recommended loading height	< 95 cm

Specifications are rounded to the nearest whole number.
 MEDIROL reserves the right to change specifications without notice.



1



2



3

Accessories

Hoverboard ViviBase

Air-sprung table with pneumatic lateral displacement



Vivera Platform S101

Light-weight table with lateral movement displacement

