

# Clubman K118

**MEDIROL**

**With 10G crash-test certification, the Clubman chair enables caregivers to safely transport patients outside and inside the ambulance without patient transfer.**

## Key features and benefits

- Anatomic patient surface for high level of comfort
- Rear telescopic lift handles
- Extendable foot-end lift handles
- Foldable side armrests
- Foldable foot support
- One four-point shoulder harness restraint
- Brakes on front or rear wheels
- Ergonomic design
- Oversized wheels maximize mobility on all types of terrain
- High visibility orange powder-coated frame
- Light-weight aluminum construction
- Pressure-washable



**250 kg**

MAXIMUM LOAD CAPACITY



**16 kg**

WEIGHT



**15 cm**

FOOT-END WHEELS DIAMETER



## Specifications

Type	<b>K118</b>
Length	<b>84 cm</b>
Width	<b>57 cm</b>
Height	<b>117 cm</b>
Weight <sup>1</sup>	<b>16 kg</b>
Wheels	
Fixed wheels diameter	<b>15 cm</b>
Swivel wheels diameter	<b>12,5 cm</b>
Maximum load capacity	<b>250 kg</b>
Recommended fastener system <sup>2</sup>	<b>F202</b>

<sup>1</sup> Chair is weighed without mattress and restraints.

<sup>2</sup> Meets dynamic crash standards EN-1789 with F202 fastening system.

Specifications are rounded to the nearest whole number.

MEDIROL reserves the right to change specifications without notice.



## Fastener

### Clubfix F202

Meets dynamic crash standards EN-1789

