

Rolman K312

MEDIROL

ROLMAN K312 is a light-weight folding chair with a stair system for easy and safe transport of patients down the stairs.

Key features and benefits

- Belted track system for movement on the stairs
- Extendable upper control handle
- Front swivel wheels
- Compact storage size
- Removable one-piece seat
- Ergonomic design
- Oversized wheels maximize mobility on all types of terrain
- High visibility orange powder-coated frame
- Light-weight aluminum construction
- Pressure-washable
- One chest restraints
- One foot restraint



180 kg

MAXIMUM LOAD CAPACITY



16 kg

WEIGHT



15 cm

BACK WHEELS DIAMETER



Specifications

| | |
|--|---------------|
| Type | K312 |
| Depth | 82 cm |
| Folded depth | 21 cm |
| Width | 55 cm |
| Height | |
| Minimum | 114 cm |
| Maximum | 155 cm |
| Weight | |
| Model without belted track system | 9 kg |
| Model with belted track system | 13 kg |
| Maximum load capacity | 180 kg |
| Recommended fastener system ¹ | F301 |

¹ Meets dynamic crash standards EN-1789 with F301 fastening system.

Specifications are rounded to the nearest whole number.

MEDIROL reserves the right to change specifications without notice.

Fastener

Unifix F301

Meets dynamic crash standards EN-1789

