

# Sanero M201

MEDIROL

**Monobloc stretcher with seven height positions and leg locking system that provides simple and ergonomic operation during height positioning and unloading from the ambulance.**

## Key features and benefits

- Steer lock system enables wheels either to swivel or to lock into straight-line position
- Seven height positions
- Automatic locking of legs when lifting
- One-person operation
- Roll-in loading with two-stage load release mechanism
- Foot-end wheel locks
- Reflective labeling
- Oversized wheels maximize mobility on all types of terrain
- Light-weight aluminum construction
- One-hand release, fold down siderails
- One-hand release, infinite positioning, lift assisted, pneumatic backrest and legrest



**250 kg**

MAXIMUM LOAD CAPACITY



**7**

HEIGHT POSITIONS



**20 cm**

WHEELS DIAMETER



## Specifications

Type	<b>M201</b>	Articulation	
Length	<b>194 cm</b>	Backrest	<b>0° - 90°</b>
Width	<b>56 cm</b>	Legrest	<b>0° - 30°</b>
Height		Headrest <sup>2</sup>	<b>± 35°</b>
Highest position	<b>89 / 100 cm</b>	Kneerest <sup>2</sup>	<b>0° - 50°</b>
Lowest position	<b>32 cm</b>	Maximum load capacity	<b>250 kg</b>
Weight <sup>1</sup>	<b>42 kg</b>	Loading height	<b>59 or 70 cm</b>
Wheels		Fastener system <sup>3</sup>	<b>F102</b>
Diameter	<b>20 cm</b>		
Width	<b>5 cm</b>		


<sup>1</sup> Stretcher is weighed without mattress and restraints. <sup>2</sup> Optional. <sup>3</sup> Meets dynamic crash standards EN-1789 with F102 fastening system. Specifications are rounded to the nearest whole number. MEDIROL reserves the right to change specifications without notice.

**EN 1789**   **EN 1865**   **CE**

### Fastener

#### Rolfix F102

Meets dynamic crash standards EN-1789



**VISIT MEDIROL.CZ TO VIEW THE ENTIRE MEDIROL ACCESSORIES**

