

Clinic Extero N114-P113

MEDIROL

Emergency patient handling system with removable stretcher options provide versatility in patient transport.

Key features and benefits

- Steer lock system enables wheels either to swivel or to lock into straight-line position
- Seven height positions
- Automatic locking of legs when lifting
- One-person operation
- Roll-in loading with two-stage load release mechanism
- Foot-end wheel locks
- Reflective labeling
- Oversized wheels maximize mobility on all types of terrain
- Extendable one-hand release lift handles
- Light-weight aluminum construction
- One-hand release, fold down siderails
- Infinite positioning, lift assisted, pneumatic backrest and legrest



250 kg

MAXIMUM LOAD CAPACITY



7

HEIGHT POSITIONS



20 cm

WHEELS DIAMETER



Stretcher specification

Type	N114
Length	195 cm
Width	56 cm
Height	12 cm
Weight ¹	20 kg
Maximum load capacity	250 kg
Handle extension	27 cm
Articulation	
Backrest	0° - 90°
Legrest	0° - 30°
Headrest ³	± 35°
Kneerest ³	0° - 50°

Transporter specification

Type	P113
Length	194 cm
Width	56 cm
Height	
Highest position	89 / 100 cm
Lowest position	32 cm
Weight	28 kg
Wheels	
Diameter	20 cm
Width	5 cm
Maximum load capacity	275 kg
Loading height	59 or 70 cm
Fastener system ²	F102


¹ Stretcher is weighed without mattress and restraints. ² Meets dynamic crash standards EN-1789 with F102 fastening system. ³ Optional. Specifications are rounded to the nearest whole number. MEDIROL reserves the right to change specifications without notice.

EN 1789 EN 1865 CE

Fastener

Rolfix F102

Meets dynamic crash standards EN-1789



VISIT [MEDIROL.CZ](https://www.medirol.cz) TO VIEW THE ENTIRE MEDIROL ACCESSORIES

