

Vivera Clinic N114-P401

MEDIROL

The unique combination of a powered transporter and detachable stretcher provides the possibility of transporting the patient even in less accessible places.

Key features and benefits

- Battery-powered hydraulic lift system
- Integrated battery with battery management system
- Automatic battery charging in fastener system
- Settable load height
- High visibility powder-coated frame
- Integrated bumper system
- One-hand release, infinite positioning, pneumatically assisted backrest and legrest
- Rubber wheels for better shock absorption
- Steer lock system enables wheels either to swivel or to lock into straight-line position
- Foot-end wheel locks



300 kg

MAXIMUM LOAD CAPACITY



5 Ah

BATTERY CAPACITY



Up to 95 cm

LOADING HEIGHT



Stretcher specification

Type	N114
Length	195 cm
Width	56 cm
Height	12 cm
Weight ³	20 kg
Maximum load capacity	250 kg
Handle extension	27 cm
Head extension ⁴	+ 29 cm
Extended patient surface ⁴	94 cm
Articulation	
Backrest	0° - 90°
Legrest	0° - 30°
Headrest ⁴	± 35°
Kneerest ⁴	0° - 50°

Transporter specification

Type	P401
Length	200 cm
Width	60 cm
Height	
Highest position	142 cm
Lowest position	47 cm
Weight ¹	71 kg
Wheels	
Diameter	15 cm
Width	5 cm
Battery capacity	5 Ah
Maximum load capacity	300 kg
Loading height	< 95 cm
Fastener system ²	F401

¹ Transporter is weighed with battery pack. ² Meets dynamic crash standards EN-1789 with F401 fastening system. ³ Stretcher is weighed without mattress and restraints. ⁴ Optional. Specifications are rounded to the nearest whole number. MEDIROL reserves the right to change specifications without notice.

IP 66 EN 1789 EN 1865 CE

Vivera Fastener F401



VISIT [MEDIROL.CZ](https://www.medirol.cz) TO VIEW THE ENTIRE MEDIROL ACCESSORIES

